

Bath Rails

MERTON MKII



A chromed or epoxy coated rail similar to the Swedish Rail but with a different shaped handle that also assists with standing up or sitting down in the bath. Length of parallel handle 175mm (7"). Length of right angled handle 250mm (10"). Height of handle 550 to 750mm (21 1/2 to 29 1/2"). Height to rim 300 to 600mm (11 3/4 to 23 1/2"). Other specifications as for AA1910 and AA1911. Weight 3.2kg.

Max user weight: 23 stone.

SWEDISH RAIL CHROME OR WHITE



This popular rail provides a firm handle just above the side of the bath, to assist with getting in and out. The handle is angled inwards to be directly over the bath rim. Available in chromed steel or epoxy coated and adjustable in height using a pin clip system.

The clamping arms can be adjusted in height and opening width. The bath rim is protected on both sides by non-slip rubber pads. The base plate of the rail must be screwed to the floor. Height of handle 547 to 747mm (21 1/2 to 29 1/2").

Height to rim 350 to 590mm (13 3/4 to 23 1/4"). Length of handle 175mm (7"). Base plate 25 x 265mm (2 x 10 1/2"). Weight 2.6kg.

BATHTUB GRAB RAIL



Providing a firm handle at the bath side, this rail helps give confidence and stability when getting in and out of the bath.

The coated steel rail is fixed to the side of the bath by a stainless steel clamp, cushioned on the inside by rubber pads.

Length 457mm (18"). Height above bath side 349mm (13 3/4"). Clamp width adjustment 83 to 140mm (3 1/4 to 5 1/2"). Weight 2.4kg

TAP RAIL



This luxury tap rail has contoured white plastic mouldings to enhance grip on the rail when raising and lowering. May be folded up when not in use. Max diameter of taps 80mm (3 1/8"). Distance of bar from taps 325mm (12 3/4"). Width of handle bar 675mm (26 1/2"). Weight 1.7kg

